


Product Datasheet

CD14 antibody (orb434888)

Description	Mouse monoclonal antibody to CD14
Species/Host	Mouse
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Chicken CD14/hulg fusion protein
Target	CD14
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Phosphate buffered saline, 0.1% Sodium Azide (NaN ₃)
Concentration	1.0 mg/ml
Storage	Keep aliquots at 2-8°C for short-term use (up to 4 weeks) and store the remaining aliquots at -20°C
Note	For research use only
Application notes	This product is suitable for use in indirect ELISA applications.
Isotype	IgG
Clonality	Monoclonal
Purity	Purified IgG prepared by affinity chromatography on Protein G from tissue culture supernatant
Expiration Date	12 months from date of receipt.



Indirect ELISA
using Mouse
anti-Chicken

...