

Product Datasheet

JAM-C antibody (orb434891)

Description	Rat monoclonal antibody to JAM-C
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IP
Immunogen	Recombinant soluble JAM-C.
Target	JAM-C
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.