

Product Datasheet

CD42a antibody (orb434913)

Description	Mouse monoclonal antibody to CD42a
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	Peripheral blood mononuclear cells.
Target	CD42a
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.