

Product Datasheet

CD3 antibody (orb434917)

Description	Mouse monoclonal antibody to CD3
Species/Host	Mouse
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FA, FACS, IHC-Fr, IP, WB
Immunogen	Porcine PBMCs
Target	CD3
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Application notes	EPSILON
Isotype	IgG1
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.