

Product Datasheet

CD14 antibody (orb434940)

Description	Mouse monoclonal antibody to CD14
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP
Immunogen	Human Thymocytes followed by peripheral blood mononuclear cells.
Target	CD14
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.