

Product Datasheet

CD335 antibody (APC) (orb434943)

Description	Mouse monoclonal antibody against CD335 conjugated to APC
Species/Host	Mouse
Reactivity	Porcine
Conjugation	APC
Tested Applications	FC
Immunogen	Fusion protein consisting of the extracellular region of porcine CD335
Target	CD335
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Allophycocyanin (APC) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.