

Product Datasheet

CD52 antibody (orb434950)

Description	Mouse monoclonal antibody to CD52
Species/Host	Mouse
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	Porcine thymocytes
Target	CD52
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Dilution Range	FACS: 1:25-200
Expiration Date	12 months from date of receipt.