

Product Datasheet

MCP-1 antibody (orb434952)

Description	Mouse monoclonal antibody to MCP-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, IHC-P, IP, WB
Immunogen	Recombinant human MCP-1.
Target	MCP-1
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.