

Product Datasheet

CD16-2 antibody (FITC) (orb434959)

Description	Rat monoclonal antibody against CD16-2 conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Fusion Protein consisting of mouse CD16-2 and rat CD4
Target	CD16-2
Preservatives	0.1% Sodium Azide (NaN ₃)1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.