

Product Datasheet

mCHERRY antibody (orb434975)

Description	Mouse monoclonal antibody to mCHERRY
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant mCherry expressed in E. coli
Target	mCHERRY
Preservatives	5% Glycerol0.1% Bovine Serum Albumin0.05mM EDTA0.03% Proclin 300. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.