

Product Datasheet

POLYUBIQUITIN (K63-LINKAGE SPECIFIC) antibody (orb434986)

Description	Mouse monoclonal antibody to POLYUBIQUITIN (K63-LINKAGE SPECIFIC)
Species/Host	Mouse
Reactivity	Aves
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, WB
Immunogen	Synthetic peptide containing the ubiquitin K63-isopeptide motif
Target	POLYUBIQUITIN (K63-LINKAGE SPECIFIC)
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.