

Product Datasheet

CD11c antibody (orb435002)

Description	Mouse monoclonal antibody to CD11c
Species/Host	Mouse
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF
Immunogen	Cells from multiple species
Target	CD11c
Preservatives	0.09% Sodium Azide (NaN ₃). Serum free tissue culture medium containing proprietary protein free supplement
Form/Appearance	Concentrated tissue culture supernatant - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Con S/N
Expiration Date	12 months from date of receipt.