

Product Datasheet

DINITROPHENOL antibody (orb435030)

Description	Rat monoclonal antibody to DINITROPHENOL
Species/Host	Rat
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Nippostrongylus brasiliensis-DNP.
Target	DINITROPHENOL
Preservatives	0.1% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgE - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgE
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.