

Product Datasheet

CD49f antibody (orb435043)

Description	Rat monoclonal antibody to CD49f
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr, IP
Immunogen	BALB/c mouse mammary tumor cells
Target	CD49f
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.