

Product Datasheet

CD11b antibody (RPE) (orb435059)

Description	Rat monoclonal antibody against CD11b conjugated to RPE
Species/Host	Rat
Reactivity	Mouse
Conjugation	RPE
Tested Applications	FC
Immunogen	Thioglycollate-elicited peritoneal macrophages (TPM)
Target	CD11b
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.