

Product Datasheet

HLA DR antibody (RPE) (orb435068)

Description	Rat monoclonal antibody against HLA DR conjugated to RPE
Species/Host	Rat
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	EHR-B Ramos cells.
Target	HLA DR
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.