

Product Datasheet

CD4 antibody (orb435105)

Description	Mouse monoclonal antibody to CD4
Species/Host	Mouse
Reactivity	Guinea pig
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr
Immunogen	Guinea pig peritoneal T-cells.
Target	CD4
Preservatives	0.1% Sodium Azide (NaN ₃) 0.7% Bovine Serum Albumin
Form/Appearance	Concentrated tissue Culture Supernatant - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Con S/N
Expiration Date	12 months from date of receipt.