

Product Datasheet

CD11b antibody (orb435111)

Description	Rat monoclonal antibody to CD11b
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P, IP
Immunogen	T cell enriched splenocytes from B10 mice.
Target	CD11b
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.