

## Product Datasheet

### CD51/CD61 antibody (orb435134)

|                            |                                               |
|----------------------------|-----------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD51/CD61        |
| <b>Species/Host</b>        | Mouse                                         |
| <b>Reactivity</b>          | Human                                         |
| <b>Conjugation</b>         | Unconjugated                                  |
| <b>Tested Applications</b> | FC, IH-Fr, IP                                 |
| <b>Immunogen</b>           | Osteoclasts from osteoclastomas.              |
| <b>Target</b>              | CD51/CD61                                     |
| <b>Preservatives</b>       | 0.09% Sodium Azide. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                         |
| <b>Concentration</b>       | 1.0 mg/ml                                     |
| <b>Storage</b>             | +4°C, -20°C if preferred                      |
| <b>Note</b>                | For research use only                         |
| <b>Isotype</b>             | IgG1                                          |
| <b>Clonality</b>           | Monoclonal                                    |
| <b>Purity</b>              | Purified                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.               |