

Product Datasheet

Goat Anti-Rabbit IgG Secondary Antibody HRP Conjugated (orb43514)

Description	Goat Anti-Rabbit IgG Secondary Antibody HRP Conjugated
Species/Host	Goat
Conjugation	HRP
Storage	Store at -20°C/1 year
Note	For research use only
Application notes	IHC/ICC: 1:500 - 1:5,000 ELISA: 1:5,000 - 1:50,000 Western blot: 1:5,000 - 1:50,000
Clonality	Polyclonal
Dilution Range	IHC/ICC: 1:500 - 1:5,000 ELISA: 1:5,000 - 1:50,000 Western blot: 1:5,000 - 1:50,000
Expiration Date	12 months from date of receipt.