

Product Datasheet

Ly-6B.2 ALLOANTIGEN antibody (FITC) (orb435141)

Description	Rat monoclonal antibody against Ly-6B.2 ALLOANTIGEN conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Cultured bone marrow cells
Target	Ly-6B.2 ALLOANTIGEN
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.