

Product Datasheet

Ly-6B.2 ALLOANTIGEN antibody (orb435144)

Description	Rat monoclonal antibody to Ly-6B.2 ALLOANTIGEN
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Cultured bone marrow cells
Target	Ly-6B.2 ALLOANTIGEN
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.