

Product Datasheet

CD3 antibody (orb435148)

Description	Mouse monoclonal antibody to CD3
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, IP
Immunogen	F344 rat T cells stimulated with PMA (TPA) and calcium ionophore
Target	CD3
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgM - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.