

Product Datasheet

CD54 antibody (FITC) (orb435156)

Description	Mouse monoclonal antibody against CD54 conjugated to FITC
Species/Host	Mouse
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Rat Ax cells (a HEV derived cell line).
Target	CD54
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.