

Product Datasheet

TUBULIN ALPHA antibody (HRP) (orb435163)

Description	Rat monoclonal antibody against TUBULIN ALPHA conjugated to HRP
Species/Host	Rat
Reactivity	Yeast
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Yeast tubulin.
Target	TUBULIN ALPHA
Preservatives	0.01% ThiomersalHRP Stabiliser (BUF052A). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.