

Product Datasheet

TUBULIN ALPHA antibody (orb435164)

Description	Rat monoclonal antibody to TUBULIN ALPHA
Species/Host	Rat
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	Yeast tubulin.
Target	TUBULIN ALPHA
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.