


Product Datasheet

CD34 antibody (orb43517)

Description	Goat polyclonal to CD34. This antibody is expected to recognize reported isoform CD34-T .
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	CD34-T
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	1-3 µg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize reported isoform CD34-T (NP_001020280.1) only. ELISA: antibody detection limit dilution 1:128000. WB: Approx 85+95kDa bands observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 35.0kDa according to NP_001764.1). The observed MW correspond to glycosylated forms. Recommended concentration: 1-3 µg/ml.
Protein Sequence	SPTGERLELEP
Clonality	Polyclonal
MW	35
Entrez	947
Dilution Range	ELISA: 1:128000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of staining of Per...



Flow cytometric analysis of Jurkat lysat...