

## Product Datasheet

### T LYMPHOCYTES antibody (FITC) (orb435176)

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody against T LYMPHOCYTES conjugated to FITC             |
| <b>Species/Host</b>        | Mouse                                                                          |
| <b>Reactivity</b>          | Rabbit                                                                         |
| <b>Conjugation</b>         | FITC                                                                           |
| <b>Tested Applications</b> | FC                                                                             |
| <b>Immunogen</b>           | Rabbit Thymocytes                                                              |
| <b>Target</b>              | T LYMPHOCYTES                                                                  |
| <b>Preservatives</b>       | 0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline          |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid |
| <b>Concentration</b>       | 0.1 mg/ml                                                                      |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                       |
| <b>Note</b>                | For research use only                                                          |
| <b>Isotype</b>             | IgG1                                                                           |
| <b>Clonality</b>           | Monoclonal                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |