

Product Datasheet

T LYMPHOCYTES antibody (orb435177)

Description	Mouse monoclonal antibody to T LYMPHOCYTES
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, IP
Immunogen	Rabbit Thymocytes
Target	T LYMPHOCYTES
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.