

Product Datasheet

CD11b antibody (orb435178)

Description	Mouse monoclonal antibody to CD11b
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IP
Immunogen	Rabbit adherent blood leucocytes.
Target	CD11b
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.