

Product Datasheet

CD45 antibody (orb435180)

Description	Mouse monoclonal antibody to CD45
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr
Immunogen	Glycoproteins isolated from the T cell line, RL-5.
Target	CD45
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.