

Product Datasheet

CD326 antibody (orb435207)

Description	Mouse monoclonal antibody to CD326
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, IHC-Fr, IHC-P
Immunogen	HT29 carcinoma cell line.
Target	CD326
Preservatives	0.09% Sodium Azide0.5%Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.