

Product Datasheet

CD326 antibody (orb435208)

Description	Mouse monoclonal antibody to CD326
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	HT29 carcinoma cell line.
Target	CD326
Preservatives	0.09% Sodium Azide0.5% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid, ready to use.
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Ready To Use
Expiration Date	12 months from date of receipt.