


Product Datasheet

NDUFA7 antibody (orb43521)

Description	Goat polyclonal to NDUFA7.
Species/Host	Goat
Reactivity	Bovine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NDUFA7 (aa27-38)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 14kDa band observed in Human Heart lysates (calculated MW of 12.6kDa according to NP_004992.2). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	RYQEISKRTQPP
Clonality	Polyclonal
MW	12.6
Entrez	4701
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...