

Product Datasheet

MACROPHAGES/MONOCYTES antibody (orb435216)

Description	Mouse monoclonal antibody to MACROPHAGES/MONOCYTES
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P
Immunogen	Human monocytes.
Target	MACROPHAGES/MONOCYTES/GRANULOCYTES
Preservatives	0.09% Sodium Azide. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.