

Product Datasheet

MACROPHAGES/MONOCYTES antibody (RPE) (orb435218)

Description	Mouse monoclonal antibody against MACROPHAGES/MONOCYTES conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Human monocytes.
Target	MACROPHAGES/MONOCYTES/GRANULOCYTES
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.