

Product Datasheet

CD45 antibody (orb435224)

Description	Mouse monoclonal antibody to CD45
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IHC-P, IP
Immunogen	Human T lymphocytes.
Target	CD45
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Low Endotoxin
Expiration Date	12 months from date of receipt.