

Product Datasheet

CD45 antibody (PerCP) (orb435229)

Description	Mouse monoclonal antibody against CD45 conjugated to PerCP
Species/Host	Mouse
Reactivity	Human
Conjugation	PerCP
Tested Applications	FC
Immunogen	Human T lymphocytes.
Target	CD45
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Peridinin Chlorophyll protein (PerCP) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.