

Product Datasheet

CD62E/CD62P antibody (orb435236)

Description	Mouse monoclonal antibody to CD62E/CD62P
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IH-Fr, IP, WB
Immunogen	Human E-Selectin (ELAM-1).
Target	CD62E/CD62P
Preservatives	0.09% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.