

Product Datasheet

CYTOKERATIN 14 antibody (FITC) (orb435250)

Description	Mouse monoclonal antibody against CYTOKERATIN 14 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FACS, ICC, IF
Immunogen	Last 15 C-terminal residues of human cytokeratin 14 conjugated to thyroglobulin.
Target	CYTOKERATIN 14
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.