

Product Datasheet

TNF ALPHA antibody (orb435262)

Description	Mouse monoclonal antibody to TNF ALPHA
Species/Host	Mouse
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr, WB
Immunogen	Recombinant ovine TNF alpha.
Target	TNF ALPHA
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.