

Product Datasheet

CD90 antibody (APC) (orb435269)

Description	Mouse monoclonal antibody against CD90 conjugated to APC
Species/Host	Mouse
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Purified human brain Thy-1.
Target	CD90
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Allophycocyanin (APC) - lyophilized
Concentration	0.1mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.