

## Product Datasheet

### CD90 antibody (orb435275)

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD90       |
| <b>Species/Host</b>        | Mouse                                   |
| <b>Reactivity</b>          | Human                                   |
| <b>Conjugation</b>         | Unconjugated                            |
| <b>Tested Applications</b> | FACS, ICC, IF, IHC-Fr, IP               |
| <b>Immunogen</b>           | Purified human brain Thy-1.             |
| <b>Target</b>              | CD90                                    |
| <b>Preservatives</b>       | None present. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                   |
| <b>Concentration</b>       | 1.0 mg/ml                               |
| <b>Storage</b>             | -20°C only (ship +4°C)                  |
| <b>Note</b>                | For research use only                   |
| <b>Isotype</b>             | IgG1                                    |
| <b>Clonality</b>           | Monoclonal                              |
| <b>Purity</b>              | Preservative Free                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.         |