

## Product Datasheet

### CD47 antibody (orb435277)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD47                            |
| <b>Species/Host</b>        | Mouse                                                        |
| <b>Reactivity</b>          | Human                                                        |
| <b>Conjugation</b>         | Unconjugated                                                 |
| <b>Tested Applications</b> | FA, FACS, IHC-Fr, IP, WB                                     |
| <b>Immunogen</b>           | Human erythrocytes                                           |
| <b>Target</b>              | CD47                                                         |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> ). TRIS buffered glycine |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                        |
| <b>Concentration</b>       | 1 mg/ml                                                      |
| <b>Storage</b>             | +4°C, -20°C if preferred                                     |
| <b>Note</b>                | For research use only                                        |
| <b>Isotype</b>             | IgG2b                                                        |
| <b>Clonality</b>           | Monoclonal                                                   |
| <b>Purity</b>              | Purified                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |