

Product Datasheet

CD55 antibody (orb435281)

Description	Mouse monoclonal antibody to CD55
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FACS, IP, WB
Immunogen	Human fibroblast cell line.
Target	CD55
Preservatives	0.1% Sodium Azide (NaN ₃) ≤ 100mM Glycine. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.