

Product Datasheet

HLA B7 antibody (RPE) (orb435323)

Description	Mouse monoclonal antibody against HLA B7 conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Papain solubilized HLA-B7 antigen
Target	HLA B7
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.