

Product Datasheet

TdT antibody (FITC) (orb435326)

Description	Mouse monoclonal antibody against TdT conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Target	TdT
Preservatives	0.1% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	50.0 ug/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Cocktail
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.