

Product Datasheet

TETANUS TOXIN antibody (orb435363)

Description	Goat polyclonal antibody to TETANUS TOXIN
Species/Host	Goat
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Purified tetanus toxin.
Target	TETANUS TOXIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.