

Product Datasheet

ESCHERICHIA COLI antibody (Biotin) (orb435367)

Description	Goat polyclonal antibody against ESCHERICHIA COLI conjugated to Biotin
Species/Host	Goat
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	A mixture of Escherichia coli serotypes.
Target	ESCHERICHIA COLI
Preservatives	0.1% Sodium azide. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.