

Product Datasheet

PTEN TUMOR SUPPRESSOR PROTEIN antibody (orb435373)

Description	Mouse monoclonal antibody to PTEN TUMOR SUPPRESSOR PROTEIN
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IP, WB
Immunogen	A 10 amino acid peptide corresponding to a region towards to the C-terminus of native PTEN
Target	PTEN TUMOR SUPPRESSOR PROTEIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.